

QXCI Channels of Information

Dr W Nelson

Many have asked what is measured by the QXCI medical device. This incredible advancement in medicine measures subtle electrical factors of the body. The history of alternative medicine has been dominated by old style resistance measuring devices. These antique devices depended on the operator to apply the probe to the acupuncture point. By subtly changing the speed of probe delivery (not the pressure), the operator could influence the outcome. This prevented true results. The other complication was the limits of just resistance testing. The body is much more complicated than just resistance. The electrical factors of the body include voltage, amperage, capacitance, inductance, frequency, and many others.

---Now don't be afraid. For over fifteen years I've been lecturing to audiences of alternative medicine, and here is where all of a sudden the apprehension and fear start. The mention of voltage, amperage etc seem to trigger fear and dread. For all of a sudden the world of electroacupuncture which before was secure and safe is now dramatically incomplete. Most of these doctors in America and Germany have spent a life time mentally invested into one channel resistance measures of energetic medicine. Now they are faced with the discovery that a new technology has so surpassed their familiarity, that at first it is intimidating. Small minds see new technology as a threat where great spirits see new technology as an opportunity.

So the QXCI represents a new technology. This article is a very very brief summary of this technology. Small minds will probably quit reading by now. Let me caution you, if you read further your life as a doctor will change and your mind might have to expand.---

The basics of bio-electrical medicine lies in the voltage, amperage, and resistance. The only thing that can truly be measured in electricity is the voltage and amperage. Everything else is a mathematical variation of voltage and amperage. So the QXCI device starts by measuring multiple channels of this information. Multiple channels are needed so that we can see variations in the electrical potential and flow of the body total. So we measure the extremities and the head. This direct measurement is of four channels of voltage, four of amperage, and four of resistance. These come into the computer for amplification and the computer acts as a frequency counter for receiving this data. This makes an active twelve channel measuring device.

Other calculations are made mathematically. Resistance can not be measured directly. It must be calculated from an amperage or voltage meter. Same as that any capacitance meter can not directly measure capacitance it must do a mathematical adjustment of voltage and amperage flow. These calculations are referred to as virtual or mathematical measures. So even the one channel resistance devices of the antique electroacupuncture devices such as listen (best), voll, phasix, bicom, mora and a host of others are virtual measures.

The QXCI device measures over forty virtual dimensions. The variations of amperage and voltage flow give us a way to measure capacitance and inductance. These are the reflections of static and magnetic effects of bio-electricity. Variations in amperage and voltage allow us to measure frequency. You can not measure a frequency with a

resistance meter. So frequency was left out of the one channel resistance systems of old.

Changes in the voltage, amperage, and resistance together make up the reactance and susceptance of the electrical system. This is the most important measure of medication testing. Resistance alone is helpful but when coupled with voltage and amperage the medication testing is greatly improved. True complete electrical reactance is the best measure of biological reactance. But multiple channels are needed to measure total reactivity.

So what is measured by the QXCI, well here are the categories measured:

1. Voltage lower left quad	33. Susceptance 1	Inverse of total
2. Voltage lower right quad	34. Susceptance 2	“ ”
3. Voltage upper left quad	35. Susceptance 3	“ ”
4. Voltage upper right quad	36. Susceptance 4	“ ”
5. Resistance lower left quad	37. Capacitance 1	Changes in amps
6. Resistance lower right quad	38. Capacitance 2	“
7. Resistance upper left quad	39. Capacitance 3	“
8. Resistance upper right quad	40. Capacitance 4	“
9. Amperage head pole 1 Brain wave	41. Capacitance 5	“
10. Amperage head pole 2 Brain wave	42. Capacitance 6	“
11. Amperage head pole 3 Brain wave	43. Capacitance 7	“
12. Amperage head pole 4 Brain wave	44. Capacitance 8	“
13. Amperage head pole 5 Brain wave	45. Inductance 1	Changes in voltage
14. Amperage head pole 6 Brain wave	46. Inductance 2	“
15. Amperage head pole 7 Brain wave	47. Inductance 3	“
16. Amperage head pole 8 Brain wave	48. Inductance 4	“
17. Frequency 1 Lower left quad	49. Inductance 5	“
18. Frequency 2 Lower right quad	50. Inductance 6	“
19. Frequency 3 Upper right quad	51. Inductance 7	“
20. Frequency 4 Upper left quad	52. Inductance 8	“
21. Frequency 5 Brain wave	53. Hydration	
22. Frequency 6 Brain wave	54. Oxidation	
23. Frequency 7 Brain wave	55. Proton an Electron pressure	
24. Frequency 8 Brain wave		
25. Reactance 1 Changes in volts, amps,		
26. Reactance 2 “ ”		
27. Reactance 3 “ ”		
28. Reactance 4 “ ”		
29. Reactance 5 “ ”		
30. Reactance 6 “ ”		
31. Reactance 7 “ ”		
32. Reactance 8 “ ”		

From these calculations several other ramifications can be surmised. The bioresonance of the system can be measured. The reaction of the body to nosodes, isodes, allersodes, sarcodes, classic homeopathics, herbals, minerals, amino acids, enzymes, hormones, bacteria, fungus, parasites, nerves, vertribrae, muscles, and thousands of other compounds. The resonance of the system is measured by determining the evoked potential or reaction of the body to these compounds.

Evoked potential is the electrical reaction of a patient to a stimulus. By analyzing the patients evoked potential or change in the above 52 measures we can see the patterns of the energetic body of the patient.

The speed of this reaction is rather quick. The body electric is drawn towards nutrition and repelled from toxins. This is discussed in detail in the Promorpheus.